POWERED
BY **Bien Air**⁺
Dental

OsseoCare™

FRA	Instructions d'utilisation
DEU	Bedienungsanleitung
ENG	Instructions for use
ESP	Instrucciones de uso
ITA	Istruzioni per l'uso
RUS	Инструкция по эксплуатации
JPN	取扱説明書
PRT	Instruções de utilização

CE
0120

Set OsseoCare™

Set OsseoCare™ REF 1700472-001

1x 1x 1x 1x 1x 1x 1x

Set OsseoCare™ REF 1700473-001

1x 1x 1x 1x 1x 1x KM 1x 1x

Set OsseoCare™ US REF 1700474-001

1x 1x 1x 1x 1x 1x 1x


Set OsseoCare™ US REF 1700475-001

1x 1x 1x 1x 1x 1x KM 1x 1x

- * - REF 1300065-001 Cable system Switzerland
- REF 1300066-001 Cable system Europe
- REF 1300067-001 Cable system USA/Asia
- REF 1306025-001 Cable system UK
- REF 1306026-001 Cable system AU

Options

Summary

Starting display

Nobel Biocare
OsseoCare™
Version 1.0

Available values

MAIN MENU	Steps	Ratio	Speed in rpm	Torque in Ncm	Irrigation in ml/min
Implantology	Drill 1	128:1	100 - 40'000 rpm with a CA 1 : 1	0.48 - 4.8 Ncm with a CA 1 : 1	30 ml/min 20%
Endodontics*	Drill 2	64:1			60 ml/min 40%
Surgery*	Tapping	30:1	Depends on the CA	Depends on the CA	90 ml/min 60%
	Tap unscrewing	27:1			120 ml/min 80%
	Implant screwing	20:1			150 ml/min 100%
	Unscrewing	16:1			
		10:1			
		1:1			
		1:2			
		1:5			



Implantology	Open pulp chamber	128:1	100 - 40'000 rpm with a CA 1 : 1	0.48 - 4.8 Ncm with a CA 1 : 1	30 ml/min 20%
Endodontics*	endo file 1	64:1			60 ml/min 40%
Surgery*	endo file 2	30:1	Depends on the CA	Depends on the CA	90 ml/min 60%
	endo file 3	27:1			120 ml/min 80%
	endo file 4	20:1			150 ml/min 100%
	endo file 5	16:1			
	endo file 6	10:1			
	endo file 7	1:1			
	endo file 8	1:2			
	endo file 9	1:5			

Implantology	Procedure 1	128:1	100 - 40'000 rpm with a CA 1 : 1	0.48 - 4.8 Ncm with a CA 1 : 1	30 ml/min 20%
Endodontics*	Procedure 2	64:1			60 ml/min 40%
Surgery*	Procedure 3	30:1	Depends on the CA	Depends on the CA	90 ml/min 60%
	Procedure 4	27:1			120 ml/min 80%
		20:1			150 ml/min 100%
		16:1			
		10:1			
		1:1			
		1:2			
		1:5			

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1 Symbols	Page
Meaning of symbols	34
2 Description	
Identification	35
Intended use	35
Environment	35
Environmental protection and information for disposal of the instrument	35
3 Set supplied	36
4 Options	36
5 Technical description	
Technical data	37
Electromagnetic compatibility	38-39
6 Installation	40-41
7 Description of keys and elements	42
8 Operation	
Description of functions	43
Start-up	44
Pre-settings (SETUP)	44-45
Description of programs	46
- Implantology	
- Endodontics*	
- Surgery*	
9 List of errors / Troubleshooting	47
10 Default values	
Implantology	48+129
Endodontics*	48+130
Surgery*	48+130
11 Maintenance	
Overhaul	48
Information	48
Cleaning-disinfection	48
Important	48
12 Generalities and guarantee	
General information	48
Terms of guarantee	48

1 Meaning of symbols

	CE Marking with number of the notified body.
	Main switch ON - The instrument is switched on. OFF - The instrument is switched off.
	Fuse Ø 5 x 20 mm.
	Alternating current.
	Device of type B.
	CAUTION ! Dangerous voltage.
	Element sensitive to electrostatic discharges.
	Warning.
	CAUTION! Refer to the accompanying documents.
	Danger of pinching. Do not put your fingers in rotating parts.
	Machine washable.
	Variability in steps.
	Symbol for "Water-cooling".
	Symbol for "peristaltic pump".
	Recyclable materials.
	Recyclable electrical and electronic materials.
	Sterilisable in autoclave up to the specified temperature.
	Intermittent operating mode.
	Manufacturer.
	Light.
	Expiration date.
	Product containing phthalates.
	Disposable product.

2 Description

Identification Electronically controlled tabletop device allowing operation of an MX-i LED micromotor with variable speed control by a pedal.
A peristaltic pump conveys the physiological liquid via a disposable irrigation line without being contaminated.
The device's LCD display indicates the stage of implant fitting, the instrument's ratio, the bur speed, torque value and irrigation flow setting.

Intended use Product intended for professional use only.
Use in dentistry for implantology, dental surgery* and endodontics* work.
Any use other than that for which this product is intended is unauthorised and may be dangerous.
The medical device meets all the current legal requirements.
The device is not authorised for use in an explosive atmosphere (anaesthetic gas).

Environment	Working	Temperature:	+10°C (50°F) to +25°C (77°F)
		Relative humidity:	30% to 80%, including condensation
		Atmospheric pressure:	700 hPa to 1060 hPa
	Transport and storage	Environmental conditions for a period of max. 15 weeks	
		Temperature:	-25°C (-13°F) to +70°C (158°F)
		Relative humidity:	10% to 100%, including condensation
		Atmospheric pressure:	500 hPa to 1060 hPa

Environmental protection and information for disposal of the instrument

The disposal and/or recycling of materials must be performed in accordance with the legislation in force.

This device and its accessories must be recycled.
Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards. The user must return the device to its dealer or establish direct contact with an approved body for treatment and recovery of this type of equipment (European Directive 2002/96/EC).

3 Set supplied

	Set OsseoCare™ REF 1700472-001	Set OsseoCare™ REF 1700473-001	Set OsseoCare™ US REF 1700474-001	Set OsseoCare™ US REF 1700475-001
OsseoCare™ Control REF 1600871-001	1x	1x		
OsseoCare™ Control REF 1600872-001			1x	1x
CA 20:1 L Micro-Series (light) REF 1600873-001	1x		1x	
CA 20:1 L Micro-Series KM (light) REF 1600874-001		1x		1x
Micromotor MX-i LED REF 1600875-001	1x	1x	1x	1x
Cable for MX-i LED Micromotor REF 1600606-001	1x	1x	1x	1x
Pack of 10 disposable sterile lines REF 1500984-010 REF 1501635-010 for KM	1x	1x	1x	1x
10 attachment collars for fastening the sterile irrig. line REF 1303711-010	1x	1x	1x	1x
Bracket for fluid bottle REF 1303393-001	1x	1x	1x	1x
3-button pedal REF 1600631-001	1x	1x	1x	1x
Cable system Switzerland, length 2.00 m REF 1300065-001	1x	1x	1x	1x
Cable system Europe, length 2.50 m REF 1300066-001	1x	1x	1x	1x
Cable system US/Asia, length 2.00 m REF 1300067-001	1x	1x	1x	1x
Cable system UK, length 2.00 m REF 1306025-001	1x	1x	1x	1x
Cable system AU, length 2.00 m REF 1306026-001	1x	1x	1x	1x
Maintenance Instructions REF 2400283	1x	1x	1x	1x
Instructions for use REF 2100259	1x	1x	1x	1x

4 Options

Contra-angle handpiece CA 20:1 L Micro-Series (light)	REF 1600873-001
Contra-angle handpiece CA 20:1 L KM Micro-Series (light)	REF 1600874-001
MX-i LED Micromotor	REF 1600875-001
Cable for MX-i LED micromotor	REF 1600606-001
3-button pedal	REF 1600631-001
Pack of 10 disposable sterile lines	REF 1500984-010
Kirschner/Meyer type detachable irrigation set for CA 20:1 L KM Micro-Series, comprising 10 rings and 10 tubes	REF 1501621-010
Kirschner/Meyer pack of 10 disposable sterile lines	REF 1501635-010
10 attachments collars for fastening the sterile irrigation line to a cable	REF 1303711-010
Bracket for fluid bottle	REF 1303393-001
Support	REF 1301575-001
Cable system Switzerland, length 2.00 m	REF 1300065-001
Cable system Europe, length 2.50 m	REF 1300066-001
Cable system US/Asia, length 2.00 m	REF 1300067-001
Cable system UK, length 2.00 m	REF 1306025-001
Cable system AU, length 2.00 m	REF 1306026-001
10x Fuse T4.0A L 250 VAC breaking capacity 40A	REF 1301560-010

5 Technical Description: Technical data

Manufacturer

Bien-Air Dental SA.

Voltage

100 – 240 VAC, 50 / 60 Hz

Fuses

2 fuses T4.0A L 250 VAC, breaking capacity 40A

Power demand

- 100 V / 300 VA
- 240 V / 300 VA

Classification

Class IIa in accordance with European Directive 93/42/EEC concerning medical devices.

Electric insulation class

I type BF, according to EN 60601-1 (apparatus protected against electric shocks).

Degree of protection

IP 40 (protection against insertion of objects larger than 1 mm).

Dimensions L x W x H

309 x 220 x 123 mm. Height with bracket 506 mm

Weight

Housing	2.7 kg	Pedal	830 g
Cable	105 g	Bracket	115 g

Memory

- Implantology mode: Storage in memory of 8 implant fitting sequences of 10 steps each.
- Endodontics* mode: Storage in memory of an endodontics sequence of 10 steps.
- Surgery* mode: Storage in memory of 4 separate programs.

Languages

French, German, English, Italian, Spanish, Portuguese, Japanese and Russian.

List of errors & Troubleshooting

Page 47

Bracket for physiological liquid flask

Stainless steel

Important

see instructions

MX-i LED Micromotor	REF 2100264
MX-i/MX-i LED micromotor cable	REF 2100163
Contra-angle CA 20:1 L, light	REF 2100263
Contra-angle CA 20:1 L KM Micro-Series, light	REF 2100263

Peristaltic pump

- Pump delivery: From 30 to 150 ml/min. (5 levels).
- Hose for pump: External Ø 5.60 mm,
Internal Ø 2.40 mm
Wall thickness 1.60 mm.

Pedal

REF 1600631-001

Dimensions (LxWxH) 250 x 205 x 54 mm
with handle: 250 x 205 x 144 mm

The pedal is waterproof (IP X8 in accordance with CEI 529).

Cables

- Length of cables:
- Pedal cable 2.90 m
- Motor cable 2.00 m

WARNING

To prevent any risk of electric shock, this device should be connected only to a power supply network provided with protective earth.

Modification of the device forbidden.

The apparatus is not adapted to the presence of an inflammable gas.

Do not attempt to open the apparatus when it is connected to the electric mains. Beware of electric shocks.

Parts applied

Micromotor MX-i LED	REF 1600875-001
Cable for MX-i LED micromotor	REF 1600606-001
CA 20:1 L Micro-Series	REF 1600873-001
CA 20:1 L KM Micro-Series	REF 1600874-001
Irrigation lines	REF 1500984-010
Irrigation lines KM	REF 1501635-010
Detachable irrigation set KM	REF 1501621-010

Accuracy of displayed values

- Speed \pm 5% (with CA Bien-Air)
- Torque \pm 5% (with CA Bien-Air)

Operating mode:

- Intermittent
- ON: 5 min.
- OFF: 40 min.

* only REF 1700472-001 / REF 1700473-001

5 Technical Description: Electromagnetic compatibility

Precautions regarding Electromagnetic Compatibility (EMC)

Electro-medical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this document.

OsseoCare™ complies with the EMC requirements according to IEC 60601-1-2. Radio transmitting equipment, cellular phones, etc. shall not be used in close proximity of the device since they could influence the performance of the device. Particular precaution is required when using strong emission sources such as High Frequency surgical equipment and similar equipment so that the HF cables are not routed on or near the device. If in doubt, please contact a qualified technician or Bien-Air Dental.

OsseoCare™ should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, OsseoCare™ should be monitored to verify normal operation in the configuration in which it will be used.

WARNING!

The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by Bien-Air Dental as replacements parts for internal components, may result in increased emissions or decreased immunity of OsseoCare™.

Guidance and manufacturer's declaration - electromagnetic emissions

OsseoCare™ is intended for use in the electromagnetic environment specified below. The customer or the user of OsseoCare™ should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	OsseoCare™ uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	OsseoCare™ is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Compliant	

Guidance and manufacturer's declaration - electromagnetic immunity

OsseoCare™ is intended for use in the electromagnetic environment specified below. The customer or the user of OsseoCare™ should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for lines no input/output	±2 kV for power supply lines ±1 kV for lines no input/output	Mains power quality should be that of a typical commercial or hospital environment.
Shock waves IEC 61000-4-5	±0.5 kV line to line ±1 kV line to line ±0.5 kV line to earth ±1 kV line to earth ±2 kV line to earth	±0.5 kV line to line ±1 kV line to line ±0.5 kV line to earth ±1 kV line to earth ±2 kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips and outages IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of OsseoCare™ requires continued operation during power mains interruptions, it is recommended that OsseoCare™ be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Essential performance: The essential performance is the maintaining of the visual lighting intensity of the LED and the maintaining of motor speed. Maximum allowed speed deviation is ± 5%.

5 Technical Description: Electromagnetic compatibility

Guidance and manufacturer's declaration - electromagnetic immunity

OsseoCare™ is intended for use in the electromagnetic environment specified below.

The customer or the user of OsseoCare™ should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of OsseoCare™, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the OsseoCare™ is used exceeds the applicable RF compliance level above, the OsseoCare™ should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the OsseoCare™.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the OsseoCare™

The OsseoCare™ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the OsseoCare™ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the OsseoCare™ as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

6 Installation


fig. 1


fig. 2


fig. 3


fig. 4

fig. 5

Installation

A. OsseoCare™ may be positioned on a table, on a trolley or another surface, but under no circumstances on the floor.
Power plug 14 is the device used for disconnection in case of problems, and it must be easily accessible at all times.

B. The fuse box may be opened with a screwdriver.
100 - 240 Vac = fuse T-4.0 A L 250 VAC REF 1301560-010

C. The equipment is powered by your line voltage (100/115/230 Vac).
Connect the power cable to the plug **fig. 1**.

D. Connect the pedal cable to the output provided on the rear panel, guiding the connector and plug by means of the index pin on the connector **fig. 2**.

! Do not lift the pedal using the connection cable.

E. Connect the micromotor cable to the motor output, guiding the connector and plug by means of the index pin on the connector **fig. 3**.

F. Align and attach the bracket to the housing provided on the rear of the console and suspend the flask or bottle **fig. 4**.

G. Check the packaging integrity, as well as the expiry date of the irrigation line.

! Warning!

Only lines supplied by Bien-Air Dental ensure trouble-free operation. These lines are sterile and for single use. Re-use may result in microbiological contamination of the patient.

AAAA/MM

6 Installation

fig. 6

H. Remove the single-use sterile irrigation line from its pouch.


fig. 7

Fitting on the spray tube

I. Connect the flexible hose of the irrigation line to the spray tube of the handpiece or contra-angle **fig. 7**.


fig. 8

Installation on the peristaltic pump

J. Install the plastic cassette in the peristaltic pump.

Check that the cassette is clipped correctly.

Close the pump lid, **fig. 8**.

If there is resistance to closing, open the lid again and check the correct positioning of the cassette.

Warning!

Do not run the pump while the lid is open.

Danger of pinching!


fig. 9

K. Perforate the cap of the physiological liquid flask with the pointed end of the irrigation line after removing the protective cap.


L. Attach the irrigation line on the motor cable using the attachment collars REF 1303711-010 **fig. 9**.

Stopping procedure

The device can be safely stopped using the main switch **16**.

7 Description of keys and elements

1 DEVICE SCREEN

2 CONTROL TO REVERSE THE ROTATION OF THE MICROMOTOR

3 "SETUP" MENU CALL-UP KEY

4 PARAMETERS BACKUP KEY

5 RETURN KEY

6 DEVICE CONTROLS

- Down key } cursor
- Up key }
- Left key (-) } adjusting
- Right key (+) }
- OK Confirmation/selection key
- In implantology/endodontics* mode: next stage

7 PERISTALTIC PUMP LID

8 MICROMOTOR CONNECTOR

9 IRRIGATION ON/OFF CONTROL BUTTON ON PEDAL

10 "PROGRAM" BUTTON ON PEDAL

- In implantology/endodontics* mode:
- Short press: next stage
- Long press: previous stage

11 BUTTON TO REVERSE THE ROTATION OF THE MICROMOTOR ON PEDAL

12 VARIABLE SPEED DRIVE ON PEDAL

13 PEDAL CONNECTOR

14 MAINS CONNECTOR (100/115/230 VAC)

15 FUSE HOLDER

16 MAIN SWITCH ON THE DEVICE

17 LABEL

18 BRACKET SUPPORT


19 MICROMOTOR

8 Operation

Description of functions

- 2 REV** The "reverse" function can be chosen directly in all the programs except endodontics*. Upon selection, a beep and the "reversal of motor rotation" icon  indicate reverse rotation. For endodontics*, settings must be performed under SETUP (Auto-forward/Auto-reverse).
- 3 SETUP** See next page.
- 4 SAVE** Stores the settings of a program: press the key until a beep is emitted, and the values that are flashing will be stored in the memory directly.
- 5 ESC** Return function. With "ESC", you can leave the current screen. In "Implantology" and "Endodontics*" mode, can also be used to return to the previous stage. If the name of the program flashes when exiting, the changes will not be taken into consideration. The changes must always be confirmed with "SAVE", otherwise they will be lost.

8 Operation

Start-up

1. System loading.....
2. LANGUAGE

English	✓
Français	
Deutsch	
Italiano	
Espanol	
Português	
Russian	
Japanese	

 Select with ▼ ▲
 Confirmation with OK
3. IMPLANT. PRG

Implant PRG1	
Implant PRG2	
Implant PRG3	✓
Implant PRG4	
Implant PRG5	
Implant PRG6	
Implant PRG7	
Implant PRG8	

 Select with ▼ ▲
 Confirmation with OK
4. <Name of selected system>
 Please check the pre-programmed values before inserting implant.
 OK: continue
 ESC: back
5. INFORMATION
 Key functions:
 ▲ move cursor UP
 ▼ move cursor DOWN
 ◀ decrem./disable
 ▶ increm./enable
 OK: continue
6. MAIN MENU

Implantology
Endodontics*
Surgery*

 Select with ▼ ▲
 Confirmation with OK
 with OK: go directly to pre-setting with no possibility of deactivating the stages
 or with SETUP, possibility of deactivating the stages with ☒


! All the pre-programmed settings are indicative and **MUST** be validated by the user.
 Memory storage takes place only when the device is first connected and is subsequently maintained.
 These parameters can then be modified in SETUP.

Pre-settings

- SETUP**
- | | |
|------------------|-----------|
| Language | English |
| Implant. PRG | Français |
| Endo system* | Deutsch |
| Ratio | Italiano |
| Light | Español |
| Footpedal | Português |
| Torque units | Russian |
| Contrast | Japanese |
| Editor | |
| System info | |
| Restore defaults | |
- Select the desired language
 ▼ ▲ and confirm with OK.
- | | |
|------------------|-----------------|
| Language | |
| Implant. PRG | |
| Endo system* | Endo steps |
| Ratio | Endo parameters |
| Light | |
| Footpedal | |
| Torque units | |
| Contrast | |
| Editor | |
| System info | |
| Restore defaults | |
- Select with ▼ ▲ then OK.
 If modifications with ◀ ▶, the screen flashes.
 save with SAVE, continue with OK.
 Without modification, continue with OK.
 ESC: back
- | | |
|------------------|-------|
| Language | |
| Implant. PRG | |
| Endo system* | |
| Ratio | 128:1 |
| Light | 64:1 |
| Footpedal | 30:1 |
| Torque units | 27:1 |
| Contrast | 20:1 |
| Editor | 16:1 |
| System info | 10:1 |
| Restore defaults | 1:1 |
- Select the desired system
 ▼ ▲ and confirm with OK.
 ESC: change
- | | |
|------------------|--|
| Language | |
| Implant. PRG | |
| Endo system* | |
| Ratio | |
| Light | |
| Footpedal | |
| Torque units | |
| Contrast | |
| Editor | |
| System info | |
| Restore defaults | |
- Select the ratio to be changed ▼ ▲ OK, then change of value with ▼ ▲ and with ◀ ▶ then save with SAVE. Continue with OK. ESC: back


8 Operation

Pre-settings

8 Operation

Description of programs

MAIN MENU	Steps	Transmission ratio	Speed in rpm	Torque in Ncm	Irrigation in ml/min
Selection wanted ▲ ▼ and confirm with OK	Each of these stages can be activated or deactivated from the SETUP menu. See also info on the last page. OK: next step ESC: previous step	Select with the cursor ■			
		Adjustable: ◀ ▶ then SAVE	Adjustable: ◀ ▶ then SAVE	Adjustable: ◀ ▶ then SAVE	Adjustable: ◀ ▶ then SAVE

Implantology	DRILL 1
Endodontics*	DRILL 2
Surgery*	TAPPING
	TAP UNSCREWING
	IMPLANT SCREWING
	UNSCREWING

Implantology	OPEN PULP CHAMBER
Endodontics*	ENDO FILE 1
Surgery*	ENDO FILE 2
	ENDO FILE 3
	ENDO FILE 4
	ENDO FILE 5
	ENDO FILE 6
	ENDO FILE 7
	ENDO FILE 8
	ENDO FILE 9

Implantology	Procedure 1
Endodontics*	Procedure 2
Surgery*	Procedure 3
	Procedure 4

128:1	100 - 40'000 rpm with a CA 1 : 1	0.48 - 4.8 Ncm with a CA 1 : 1	30 ml/min 20%
64:1			60 ml/min 40%
30:1			90 ml/min 60%
27:1	Depends on the CA selected	Depends on the CA selected	120 ml/min 80%
20:1			150 ml/min 100%
16:1			
10:1			
1:1			
1:2			
1:5			

Pre-programmed values, see pages 129-130

9 List of errors & Troubleshooting

Message		Cause of error	Action
 Release the pedal	The pedal is pressed when starting the device. The motor is blocked for more than 2 sec.	Safety	Release the pedal and press again.
 	The motor control card limits the power supplied to the motor to prevent motor overheating.	Safety	Avoid extended use.
Equipment initialisation error The following error may occur at start-up of OsseoCare™			
1. Check on the integrity of the OsseoCare™ memory			
INIT ERROR 1	The memory is corrupt! Please contact Bien-Air Dental SA. ESC: restore	The memory data check failed.	Press the ESC key to try to restore the memory. Contact Bien-Air Dental SA.
Device operating error The following errors may occur during operation of the device			
1. Loss of pedal connection			
ERROR 1	The pedal is not connected! Please check the connection. ESC: exit	The pedal is not correctly connected.	Check pedal connection. Contact Bien-Air Dental SA.
2. Peristaltic pump overheating			
ERROR 2	Irrigation pump overheating! Please wait for it to cool. ESC: exit	Peristaltic pump motor overheating	Wait until the system cools. Contact Bien-Air Dental SA.
3. Peristaltic pump general error			
ERROR 3	Irrigation pump fault! Please contact Bien-Air Dental SA. ESC: exit	Peristaltic pump electrical fault.	Contact Bien-Air Dental SA.
4. Loss of motor connection			
ERROR 4	The motor is not connected! Please check the connection. ESC: exit	Loss of motor phase fault. The motor is not correctly connected.	Check motor connection. Contact Bien-Air Dental SA.
5. Motor cable fault			
ERROR 5	Motor cable fault! Please change cable. ESC: exit	Motor power fault. The motor cable may be defective.	Check motor cable. Contact Bien-Air Dental SA.
6. Motor control overheating			
ERROR 6	System overheating! Please wait for it to cool. ESC: exit	Overheating of motor control card (electrical control of motor).	Wait until the system cools. Contact Bien-Air Dental SA.
7. System electrical fault			
GEN ERROR [Error code]	System electrical fault! Please contact Bien-Air Dental SA. ESC: exit	Communication fault with motor control card: [EC100] Motor control card power supply undervoltage: [EC101] Motor control card power supply overvoltage: [EC102] Other motor control card faults: [EC120]	Contact Bien-Air Dental SA.

10 Default values

Implantology:

Default values page 129

The table shows the default operating values for the 8 implantology systems pre-programmed in the system.

Endodontics*:

Default values page 130

The table shows the default operating values for the endodontics sequence.

11 Maintenance

Only use original Bien-Air Dental maintenance products and parts or those recommended by Bien-Air Dental. Using other products or parts may cause operational failure and/or void the guarantee.

Overhaul

Never disassemble the device. For any modification and repair, we recommend that you contact your regular supplier or Bien-Air Dental directly. Bien-Air Dental asks the user to have its dynamic instruments checked or inspected at least once a year.

Information

The technical specifications, illustrations and dimensions contained in these instructions are given only as a guide. They may not be the subject of any claim. The manufacturer reserves the right to make technical improvements to its equipment, without amending these instructions. For all additional information, please contact Bien-Air Dental SA at the address indicated on the back cover.

12 Generalities and guarantee

General information

The device must be used by a qualified person in compliance with the current legal provisions concerning industrial safety, health and accident prevention measures, and these working instructions. In accordance with such requirements, the operators:

- must only use devices that are in perfect working order; in the event of irregular functioning, excessive vibration, abnormal heating or other signs that may indicate malfunction of the device, the work must be stopped immediately; in this case, contact a repair centre approved by Bien-Air Dental;
- must ensure that the device is used only for the purpose for which it is intended, must protect themselves, their patients and third parties from any danger, and must avoid contamination through the use of the product.

Terms of guarantee

Bien-Air Dental grants the user a guarantee covering all functional defects, material or production faults.

The device is covered by this guarantee for 24 months from the date of invoicing.

In case of justified claim, Bien-Air Dental or its authorised representative will fulfil the company's obligations under this guarantee by repairing or replacing the product free of charge.

Surgery*:

Default values page 130

The table shows the default operating values for 4 different surgical sequences.

Cleaning-disinfection

- Disinfect the surfaces of the console and pedal with a clean cloth soaked in a suitable product.
- Do not exert any pressure on the screen.
- Do not immerse in disinfectant solution.
- Not designed for an ultrasonic bath.
- Use a new sterile irrigation line for each patient.

Important

For maintenance:	see instructions
MX-i LED micromotor	REF 2100264
Cable for MX-i LED micromotor	REF 2100163
Contra-angle CA 20:1 L, light	REF 2100263
Contra-angle CA 20:1 L KM Micro-Series, light	REF 2100263

Any other claims, of whatever nature, in particular in the form of a claim for damages and interest, are excluded.

Bien-Air Dental shall not be held responsible for damage or injury and the consequences thereof, resulting from:

- excessive wear and tear • improper use • non-observance of the instructions for installation, operation and maintenance
- unusual chemical, electrical or electrolytic influences • poor connections, whether of the air, water or electricity supply.

The guarantee does not cover flexible "fibre optic" type light conductors, or any parts made of synthetic materials.

The guarantee shall become null and void if the damage and its consequences are due to improper manipulation of the product, or modifications to the product carried out by persons not authorised by Bien-Air Dental. Claims under the terms of the guarantee will be considered only on presentation, together with the product, of the invoice or the consignment note, on which the date of purchase, the product reference and the Serial No. should be clearly indicated.

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